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“New techniques for diagnosis of diatom *Chaetoceros* and *Rhizosolenia*”

19-30 May 2014: Techniques for isolation of single cells, establishing and culturing species, and DNA extraction methods

- **19-21 May 2014:** Single cell isolations for establishing cultures of type species of diatoms from Thailand, maintained in L1-medium (Guillard&Hargraves, 1993) at a salinity of 31 and a temperature of 22-26°C
- **22-23 May 2014:** Learn further techniques for isolation of single diatom cells from field samples. Discussion about different culturing techniques.
- **26-27 May 2014:** Harvesting and concentrating cellular material from the cultures, to a frozen and stored for DNA analyses. This is to be continued during the following weeks.
- **28-30 May 2014:** Learn more about DNA extraction methods and precipitation of DNA using different methods

1-30 June 2014: Cleaning DNA, method for PCR-product

- **2-6 June 2014:** Learn more techniques for cleaning DNA and PCR amplification of selected genes. Meeting for discussions about selection of genes for different purposes, problems with techniques and culturing
- **9-13 June 2014:** Learning and practicing methods for single cell PCR amplification of diatoms from field samples. Sequencing of PCR products. Start analyzing data and writing paper
- **16-20 June 2014:** Further morphological and molecular studies on established cultures and single cells from field samples. Discussions about different types of phylogenetic analyses, and their advantages and drawbacks

- 23-27 June 2014: Continued discussion and practising of phylogenetic analyses. Discuss data and potential problems. Searching useful references from the Royal Library of Copenhagen
- 30 June 2014: Finalizing reference search from the Royal Library of Copenhagen and the Botanical Library

1-8 July 2014: Try to get all result for write the paper for publication

- 1-4 July 2014: Continued analyses and writing of paper
- 7 July 2014: Discussion of data and other results for publication and plans for the future
- 8 July 2014: Gathering and checking data for continued work in Thailand

DNA extraction

Equipment

1. Ice tray
2. 2x CTAB, Chloroform, Isoamyl alcohol, 96% & 70%EtOH
3. 3M sodium acetate (NaAc)
4. Temp's centrifuge (in hood)
5. Oven (60-65 °C), Vortex machine
6. Small centrifuge machine
7. Pipette (200P and 1000P) and tip
8. 1.5 ml micro-tubes = no.sample x 3
9. Parafilm, beaker and erlenmeyer flask 20 ml

Method

1. Keep CTAB (from refrigerator) in 60 °C oven until disappear small white ball but from the table it already!


2. Centrifuge samples (small centrifuge) few min and throw away all water from samples and keep tubes in rack on ice tray.

3. Run no. of samples on 1.5 ml micro-tubes (1).
4. Add **500 μ l CTAB** in all 3. add oven 10 min if want add iron ball for break chain add here 2 min
5. Add few samples in 4. (be careful! same name on each tube) and vortex
6. Keep in 60 $^{\circ}$ C oven and vortex every 10 min in 1 hr.(about 5-6 time)
7. Make 2 set 1.5 ml micro-tubes code set (no.2) and real name set (no.3)

Rinse with Chloroform

1. Make 960 μl Chloroform + 40 μl isoamyl alcohol (Isopentyl ethanol) for all samples
(calculate 24 Chloroform : Isoamyl alcohol see in table 1)

2. Add 500 μl 24:1 of Chloroform:Isoamylalcohol to every tubes and turn tubes
3. Centrifuge 10,000 G 10 min Temp. 4 $^{\circ}\text{C}$ (centrifuge in hood)

Clean iron ball: add 0.4M HCl 20 min
and add 96% EtOH and leave in oven

4. Take supernatant (upper layer) to new micro-tube no.2 (P1000 and P200)
5. Add 500 μl 24:1 of Chloroform:Isoamylalcohol to every tubes and turn tubes

6. Centrifuge **10,000 G 10 min** Temp. 4°C (centrifuge in hood)
7. Take supernatant to new micro-tube no.3

Precipitation of DNA (EtOH in refrigerator)

1. Add **1,000 μl** 96% EtOH and Add **50 μl** 3M Sodium acetate (NaAc) and shake
2. Centrifuge 10 min at **20,000G** 4°C
3. Remove carefully liquid without disturbing pellet (P1000)
4. Add **500 μl** 70% EtOH, wash gently and centrifuge **10 min at 20,000G**
5. Remove carefully liquid without disturbing pellet (P1000) and centrifuge **10 min at 20,000G**
6. Remove remaining liquid (P200)
7. Dry in rack at 60°C in oven (few minutes ONLY!) add **31 μl** MQ-water (distill water) and close with parafilm and leave at 30°C until next day.

Be careful:

1. Keep DNA samples on ice tray all time.
 2. Chloroform is carcinogen reagent please make in hood only.
-

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